

## 510(k) SUMMARY

### K191401 - PregSense™

#### Submitter

Nuvo-Group Ltd.  
Yigal Alon 94 St., Alon Tower 1  
Tel Aviv, Israel 6789155  
Phone: +972-54 234 7770  
Contact Person: Adar Shani

Date Prepared: March 26, 2020

**Name of Device:** PregSense™

**Common or Usual Name:** Home Uterine Activity Monitor

**Regulation Number:** 21 CFR 884.2730

**Regulation Name:** Home Uterine Activity Monitor

**Product Code:** LQK

**Product Code Name:** Home Uterine Activity Monitor

**Regulatory Class:** II

**Predicate Device:** Sense4Baby System Model B+ (K143114). The predicate device has been subject to a design-related recall (belt clip defects).

#### Device Description

PregSense™ is a non-invasive medical device that acquires and displays vital signs of the pregnant woman and of her fetus. It measures and processes signals picked up on the abdominal surface using sensors, electronic circuitry and processing software. Two types of sensors pick up the signals: electrocardiogram (ECG)-like sensors that capture bio-potential signals, and acoustic sensors. The bio-potential (ECG-like) sensors capture fECG (heartrate of the fetus) and mECG signals (heartrate of the pregnant woman). The acoustic sensors measure the sounds from the pregnant woman's abdomen, (PCG -phonocardiogram and fPCG -fetal PCG). Monitoring of the fetal and maternal heart rate using PregSense is limited to a five-minute session.

PregSense™ is an integrated platform that uses a signal acquisition tool to provide input to two separate software applications, one for the patient (PregSense™ ME) and one for the physician (PregSense™ MD). The sensors are incorporated in a belt (PregSense™ Belt) that is worn on the abdomen of the pregnant woman, where it acquires both biopotential and acoustic signals. The signals are processed at the cloud-server level where the inputs from the sensors are processed, merged and downloaded to the mobile devices of the pregnant woman and her health care provider. The PregSense™ ME application (for the patient) allows the pregnant woman to view the average maternal and fetal heart rate after a

five minute session has been completed and PregSense™-MD application (for the health care provider) allows the health care provider to view the complete fetal and maternal heart rate data from the five minute session online and remotely via the internet. A monitoring session can only be initiated by a health care provider.

**Indications for Use**

PregSense is a maternal-fetal monitor that non-invasively measures and displays fetal heart rate (FHR) and maternal heart rate (MHR). The PregSense acquires and displays the FHR and MHR tracings from abdominal surface electrodes that pick up the fetal heart biopotential and maternal heart biopotential signal, and from surface acoustic sensors that pick up the fetal PCG (fPCG; phonocardiogram) and the maternal PCG (mPCG; phonocardiogram) signals.

PregSense is indicated for use by pregnant women who need documentation of fetal heart rate activity, and who are in their 32nd week of gestation (or later), with a singleton pregnancy. PregSense is intended to be used for a maximum of five minutes.

The PregSense maternal-fetal monitor is intended for use in the antepartum period by healthcare professionals in health care facilities and by the patient in the patient’s home, on the order of a physician.

The PregSense is not intended for use in critical care situations or in laboring patients or those patients hospitalized for or suspected to have preterm labor.

PregSense is not intended to be used for antepartum monitoring (e.g., non-stress testing).

**Substantial Equivalence Comparison**

A table comparing the intended use and technological characteristics of the subject and predicate devices is provided below.

**Comparison Chart**

	<b>Subject Device: Nuvo’s PregSense System</b>	<b>Predicate Device: Sense4Baby System Model B+</b>
<b>510(k) Number</b>	K191401	K143114
<b>Product Code</b>	LQK	LQK, MOH, HGM
<b>Classification</b>	21 CFR 884.2730	21 CFR 884.2730
<b>Device Type</b>	Maternal-fetal monitor	Maternal-fetal monitor
<b>Intended Use/ Indications for Use</b>	PregSense is a maternal-fetal monitor that non-invasively measures and displays fetal heart rate (FHR) and maternal heart rate (MHR). The PregSense acquires and displays the FHR and MHR tracings from abdominal surface electrodes that pick up the	The Sense4Baby System Model B+ is indicated for conventional antepartum fetal monitoring applications in pregnancies greater than or equal to 24 weeks gestation. It may be used for antenatal monitoring (e.g., non-stress testing and/or uterine activity

	<b>Subject Device:</b> <b>Nuvo's PregSense System</b>	<b>Predicate Device:</b> <b>Sense4Baby System Model B+</b>
	<p>fetal heart biopotential and maternal heart biopotential signal, and from surface acoustic sensors that pick up the fetal PCG (fPCG; phonocardiogram) and the maternal PCG (mPCG; phonocardiogram) signals.</p> <p>PregSense is indicated for use by pregnant women who need documentation of fetal heart rate activity, and who are in their 32nd week of gestation (or later), with a singleton pregnancy. PregSense is intended to be used for a maximum of five minutes.</p> <p>The PregSense maternal-fetal monitor is intended for use in the antepartum period by healthcare professionals in health care facilities and by the patient in the patient's home, on the order of a physician.</p> <p>The PregSense is not intended for use in critical care situations or in laboring patients or those patients hospitalized for or suspected to have preterm labor.</p> <p>PregSense is not intended to be used for antepartum monitoring (e.g., non-stress testing).</p>	<p>monitoring) in a health care setting or home.</p> <p>It is to be used by health care professionals and patients on the order of a physician.</p> <p>Before the Sense4Baby System Model B+ is prescribed for home use, the user (patient) must be instructed/trained in proper use of the equipment.</p> <p>Home uterine activity monitoring has not been shown to prevent the onset of preterm labor nor will it prevent the occurrence of preterm birth.</p>
<b>Prescription Use</b>	Medical personnel or patients on order of physician	Medical personnel or patients on order of physician
<b>Intended Environments</b>	Health care setting or home	Health care setting or home
<b>Target Population</b>	Women who are $\geq 32$ gestational weeks with singleton pregnancy	Women who are $\geq 24$ weeks gestation
<b>Patient Interface</b>	Maternal abdomen connected to surface ECG-like bio-potential sensors	Maternal abdomen connected to sensor with Piezo-electric crystals

	<b>Subject Device:</b> <b>Nuvo's PregSense System</b>	<b>Predicate Device:</b> <b>Sense4Baby System Model B+</b>
	and acoustic sensors	
<b>Data Collected from Sensor Array</b>	Fetal heart rate, Maternal heart rate	Fetal heart rate, Maternal heart rate, Uterine activity
<b>Technology Employed</b>	Transabdominal electrocardiography signals and acoustic signals	Pulsed Doppler Ultrasound
<b>Monitoring Session</b>	5 minutes	30 minutes
<b>Information Displayed On</b>	Cloud or mobile based software applications	Web based portal

The subject and the predicate device have different indications for use statements, but have the same intended use – for antenatal monitoring in a health care setting or home environment by the patient by prescription only.

The subject device is to be utilized in women who are greater than 32 weeks of gestation, whereas the predicate device can be used in subjects who are greater the 24 weeks of gestation. Both the subject and predicate device are intended to be used for conventional antepartum monitoring by both the physician and the patient. The predicate device is indicated for use for both fetal/maternal heart rate monitoring and uterine activity monitoring, whereas the subject device is indicated for fetal/maternal heart rate monitoring only. However, as both the subject and predicate device are intended to be used for fetal/maternal heart rate monitoring, the lack of uterine activity monitoring does not raise different questions of safety and effectiveness as it pertains to the function of the device. In addition, the subject device is intended to be used for a maximum of five minutes, rather than 30 minutes as indicated for the predicate device. This is not a new intended use, as both the subject and predicate device are used to document fetal and maternal heart rate activity.

The subject and predicate devices have different technological characteristics. The FHR/MHR technology (i.e., ECG vs pulsed doppler ultrasound) of the subject and predicate device are different, as are the related hardware components associated with these different technologies. In addition, the subject device does not contain uterine activity monitoring technology. Finally, the method of data display between the subject and predicate device is different, with the subject device using a cloud-based service and software apps and the predicate device using a web based portal. However, different types of safety and effectiveness questions are not raised by these differences in technological characteristics.

### Non-Clinical Performance Data

The following non-clinical performance testing was provided to support the performance of PregSense™:

- Electronical testing
  - electromagnetic testing (IEC 60601-1-2),
  - electrical safety testing (IEC 60601-1),
  - battery safety
  - external defibrillation safety testing,
  - testing of the electrical interface and electronic parts,
  - Bluetooth functionality,
  - over-temperature protection, safety tests, and relevant use cases;
- Functionality testing to evaluate the durability and functionality of the various parts of the device, including hardware and accessories
  - PCB testing
  - Belt electronic testing
  - Short circuit protection testing
  - Mechanical functionality
  - Cable and rubber webbing testing
  - Velcro degradation
  - ECG acquisition/processing/detection/classification per 60601-2-27 and 606001-2-47
- Software verification and validation testing for the PregSense ME App, Software version: 0.40.1.3669; PregSense MD App, Software version: 0.2 (117); PregSense WSH firmware version: 0.123; and the PregSense Server, Software version 0.1.90 per the recommendations of the 2005 guidance document *Guidance for the Content of Premarket Submission for Software Contained in Medical Devices*
- Cleaning and disinfection information per the recommendations of the 2015 guidance document *Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling*
- Biocompatibility testing per the recommendations of the 2016 guidance document *Use of International Standard ISO 10993-1, "Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process" as follows:*
  - Cytotoxicity – ISO 10993-5:2009/®2014
  - Sensitization – ISO 10993-10:2010
  - Irritation – ISO 10993-10:2010
- Human Factors testing per the recommendations of the 2016 guidance document *Applying Human factors and Usability Engineering to Medical devices*

### Clinical Performance Data

Two clinical studies were performed using the belt with embedded sensors and associated software for analysis: a feasibility study and a pivotal study. The feasibility study was conducted to assess the feasibility and safety of PregSense™ system. The pivotal trial was conducted in order to demonstrate

that the PregSense™ performs appropriately for its intended use. The feasibility and pivotal trials were conducted using the device in thirty-minute sessions. However, the data used to support the safety and effectiveness of the subject device was limited to five minutes.

During the feasibility study, a total of 76 subjects participated for which a total of 510 recording sessions were executed. Female subjects aged  $\geq 18$  years and  $\leq 50$  years with a singleton pregnancy between 20–40 weeks of gestation and who were capable of signing informed consent were included in the study. No adverse events, procedure-related or device-related, were reported during the study. The results of the feasibility study showed that the overall detection percentage was  $> 70\%$  for the overall population and  $90\% (\pm 11.3)$  for those at least 32 weeks of gestation. Therefore, the pivotal study was performed in women at 32 weeks or more.

A pivotal clinical study in four healthcare settings (two in the United States and two outside the United States) to permit comparison to a “standard of care” monitoring device (CTG) was performed in 149 subjects. Because this type of “standard of care” monitoring could not be performed in a home environment, use in a healthcare setting was required. To ensure the device will perform appropriately when the belt is worn by the patient at home, a usability study was performed with lay users. Using Bland-Altman Limits of Agreement, the pivotal clinical testing for the complete 30-minute session demonstrated that for FHR, 97.08% of all differences lie between the 95% agreement limits of  $[-8.84, 8.24]$  bpm in the 149 subjects with no adverse events reported. In addition, the pivotal clinical testing demonstrated that for MHR, 95.31% (48241/50616) of all differences lie between the 95% agreement limits of  $[-5.30, 5.86]$  bpm.

As the duration of data collection of the PregSense system is limited to 5 minutes, an additional analysis of the limits of agreement using only 5 minutes of outputted data post-calibration showed narrower limits of agreement than originally reported using the entire 30-minute data collection: the lower limit of the 95% confidence interval of the lower agreement bound is  $-7.47$  bpm and the upper limit is  $8.16$  bpm. Furthermore, 96.64% of the differences fall between the limits of agreement.

Together, the completed clinical and usability testing demonstrate that the device accuracy is acceptable and that lay users can use the device successfully.

### **Conclusions**

The results of the performance testing described above demonstrate that the PregSense device is as safe and effective as the predicate device and supports a determination of substantial equivalence.